



Test Report

Determination of the retention capacity of
Erlab Filters for methylene chloride
according to the new EPA regulations.



Supervisor- Marie Fauvre

Erlab Research & Development Director

Author - Laurane Michaux

Erlab R&D Laboratory Technician

Specialist in air purification and gas analysis processes

SUMMARY

I. Introduction	5
II. Methylene chloride detection at TD/GC/MS	6
II.1 Overview of Gas Chromatography and Mass Spectrometer	6
II.2 Quantification of methylene chloride	9
II.3 Implementation for methylene chloride concentration measurements	12
III. Presentation of the results	15
III.1 Neutrodine Unisorb®	15
III.2 AS	17
III.3 Neutrodine®	19
III.4 Summary	21
IV. Conclusion	22

Table of illustrations

Figures

Figure 1: Comparison of OSHA's Methylene Chloride Exposure Limits with the New EPA Regulations	5
Figure 2: Gas chromatography - mass spectrometer.....	6
Figure 3: Gas chromatography analysis method	7
Figure 4: Mass spectrometer analysis method	8
Figure 5: Identification of methylene chloride using the NIST library	8
Figure 6 : Sac Tedlar de 1L	9
Figure 7: Quantification Method	10
Figure 8: Photo of the filtration test on Captair Smart 321 linked to the extraction	12
Figure 9: Schematic diagram of the hot-generation filtration test	13
Figure 10: STUART Griddle (Model SB460, Serial Number R000101212)	13
Figure 11: HEIDOLPH Peristaltic Pump (Pumpdrive Model 5101, No. 523-21010-00-2, Serial Number 111101371)	14
Figure 12: OHAUS "Adventurer" Scale (AX201 Model, Serial Number B709751985)	14

Tables

Table 1 : Summary of the tests on Neutrodine Unisorb® in class 1	17
Table 2 : Summary of tests on AS in class 1	18
Table 3 : Summary of the tests on Neutrodine® in class 1	20
Table 4 : Summary of values at 0.02 ppm.....	21
Table 5 : Summary of values at 0.5 ppm.....	21
Table 6 : Summary of times at 0 ppm	21

Graphics

Graph 1: Methylene chloride calibration (0.01 ppm – 0.5 ppm)	11
Graph 2: Methylene chloride calibration (full range 0.01 – 50 ppm)	11
Graph 3: Neutrodine Unisorb® Class 1 breakthrough curve as a function of evaporated mass.....	15
Graph 4: Neutrodine Unisorb® Class 1 Breakthrough Curve as a Function of Time	16
Graph 5: AS class 1 breakthrough curve as a function of evaporated mass.....	17
Graph 6: Breakthrough curve AS class 1 as a function of time	18
Graph 7: Neutrodine® breakthrough curve class 1 as a function of evaporated mass	19
Graph 8: Neutrodine® breakthrough curve class 1 as a function of time	20

I. Introduction

In November 2022, the U.S. Environmental Protection Agency (EPA) updated its regulations regarding the risk determination of methylene chloride. This revised version concludes that this chemical substance poses a harmful risk to health. The main health risks identified are: neurotoxicity due to short-term exposure, effects on the liver and cancer risk due to long-term exposure.

The new EPA regulations have significantly reduced the exposure threshold compared to the OSHA (Occupational Safety and Health Administration) standard, as shown by the Figure 1.

Rule	OSHA	New EPA
8-Hour Time Weighted average (TWA)	25 ppm	2 ppm
15-Minute Short Term Exposure Limit (STEL)	125 ppm	16 ppm
Action Level	12.5 ppm	1 ppm

Figure 1 : Comparison of OSHA's Methylene Chloride Exposure Limits with the New EPA Regulations

In accordance with the requirements of the NFX 15-211 standard for the protection of laboratory operators, a Class 1 filtration hood equipped with a main filter and a safety filter must ensure a maximum concentration downstream of the filters of less than 1% of the TLV of the molecules handled. To determine the retention capacity of Erlab filters according to 1% of the EPA's TVL, a specific analytical method was developed to ensure a detection limit of less than or equal to 0.02 ppm methylene chloride.

This method is based on the analysis of released methylene chloride concentration by Erlab filters by gas chromatography-mass spectrometry coupled to a thermodesorber (TD/GC/MS).¹ Throughout the test, a known concentration of methylene chloride is evaporated inside the fume hood for several hours. Air samples are continuously taken downstream of the Erlab filters and analyzed to determine the amount in grams of methylene chloride that can be evaporated in the filter hood down to a concentration of 0.02 ppm.

¹ TD/GC/MS = thermodesorber / gas chromatography / mass spectrometry

This report presents the results of filtration efficiency tests on a 2C configuration filtration hood equipped with Erlab filters. Three types of filters are tested simultaneously: 2C-GF4AS, 2C-Neutrodine Unisorb and 2C-Neutrodine.

I. Methylene chloride detection at TD/GC/MS

I.1 Overview of Gas Chromatography and Mass Spectrometer



Figure 2 : Gas chromatography - mass spectrometer

Gas chromatography coupled with mass spectrometer, illustrated in Figure 2, is an analytical method that separates chemicals from a sample, detecting and identifying its components based on their mass-to-charge ratio. This precise method allows for the analysis of many trace chemical components.

The R&D laboratory is equipped with this equipment and a thermodesorber to trap pollutants in a cold trap and inject them into the gas chromatograph.

The gas chromatograph is set to a specific method for separating VOCs from a gas sample
(Figure 3).




Figure 3 : Gas chromatography analysis method

For these tests, the mass spectrometer is configured in Selected Ion Full Ion (SIFI) mode. The SIFI mode includes the Selected Ion Recording (SIR) and Full Scan (FS) modes. The SIR mode is more sensitive and focuses on selected ions while the FS mode allows for full spectral scanning. The SIFI mode therefore allows the user to obtain the spectral information of the FS with the sensitivity of the SIR. The device has a spectral database, called NIST, used to identify compounds. SIFI mode becomes useful when it comes to obtaining information with higher sensitivity from existing GC/MS methods, without additional analyses.

During the Full Scan of the air sampling at the outlet of the filter column, the mass spectrometer is set to detect and quantify ions, thanks to their mass-to-charge ratio, denoted m/z , over a range of m/z 10 to 300 for 8.33 min. At the same time, the SIR mode detects the presence of the characteristic ions of methylene chloride based on the mass-to-charge ratio of 3 to 5 min (Figure 4).

Thanks to FS, we know that the retention time of methylene chloride is between 4 and 5 minutes.



To identify and quantify methylene chloride, the two main ions are used: $m/z = 49$ and $m/z = 84$ (Figure 5).



I.1 Quantification of methylene chloride

Methylene chloride is a substance that is easy to analyze by TD/GC/MS, because it is an organic compound. In addition, the existence of qualifying ions in the library facilitates the identification of the compound of interest.

To have an accurate quantification over a wide range of concentrations, we have chosen to create two ranges that will be used simultaneously according to the concentration of dichloromethane to be quantified.

To quantify concentrations below 0.5 ppm detected downstream of the filter, we performed a range from 0.01 ppm to 0.5 ppm. The concentrations in this range are:
0.01 ppm / 0.02 ppm / 0.05 ppm / 0.1 ppm / 0.5 ppm (Graph 1).

A second range to quantify the concentrations upstream of the filter has been carried out. The concentrations for this range are: 5 ppm / 10 ppm / 30 ppm / 50 ppm.

To make these ranges, we use 1L Tedlar bag (Figure 6), previously filled with 500 mL of clean air, and dedicated to the same product to avoid contamination.



Figure 6 : Sac Tedlar de 1L

We start with a “father” bag in which a volume of methylene chloride is injected into liquid form to obtain a precise concentration. From this “father” bag, a series of dilutions is made to obtain a series of points at increasing concentrations, used to establish the calibration curve. The measurements of these points are repeated three times to ensure the repeatability of the results.

After analysis of the bags, the quantification method (Figure 7) was created from the mid-range chromatogram and the response of the peak (area) as a function of concentration was obtained. Some points (having outliers) and the y-intercept have been excluded.

Method Editor - DCM

File Edit Help

Compound:

1: DCM

Name: DCM

Internal Ref: [None]

Data Source: ☒ Mass Spec ☐ GC-A ☐ GC-B

Quantify Trace: 49

Acquisition Function Number: Any

Concentration of Standards: Conc. A

Peak Location

☒ Retention Time (mins): 4.284

☐ Relative Retention Time: 0.000

Time Window (mins) \pm : 0.200

Peak Matching

Peak Selection: Multiple Ion Ratio-to Quantify Trace

REV Fit Threshold: 0

Buttons: Append, Insert, Modify, Delete, General Parameters..., Integrate Parameters..., Environmental Parameters

☐ User RF Value: 1.000000

User Peak Factor: 1.000000

Reporting Threshold: 0.000

Standard Concentration Factor: 1.000

Spectrum Multi-Ion

Qualifier Ions

	m/z	Target Ratio	\pm Tolerance(%)
1.	84.00	58.9	20.0
2.			20.0
3.			20.0
4.			20.0

Tolerance: Relative


Coelution Window (sec) \pm : 1.00

Figure 7 : Quantification Method


To plot the breakthrough curve and have an overall overview of the operation of the Erlab filter beyond 0.5 ppm, a so-called “complete” range grouping all the injected points is drawn (Graph 2). So, to determine the mass of methylene chloride in grams evaporated at 0.02 ppm, we will use the range 0.01 - 0.5 ppm to be as accurate as possible. And beyond this concentration, the full range of 0.01 – 50 ppm will be used.

Test Report: Determination of the retention capacity of Erlab Filters for methylene chloride

According to new EPA regulations



Graph 1 : Methylene chloride calibration (0.01 ppm – 0.5 ppm)



Graph 2 : Methylene chloride calibration (full range 0.01 – 50 ppm)

Thanks to these calibration lines, we can then quantify the methylene chloride detected upstream and downstream of the filtration column and taken as described below.

1.1 Implementation for methylene chloride concentration measurements

The test configuration is as follows (Figure 8): gas sampling is carried out at the outlet of the Captair Smart 321 hood equipped with a 2C filtration column connected to the extraction for the duration of the test until the filters are saturated. Depending on the type of filter, the duration of the test can vary from one to several days. The opening area used for these tests is the Oblong, which has an opening area of 0.12 m².

The filters tested are Neutrodine, Unisorb®, AS, and Neutrodine®.




Figure 8 : Photo of the filtration test on Captair Smart 321 related to the extraction

The analyzer is connected to a sampling grid that allows homogeneous aspiration of VOCs at the filter outlet.

Figure 9 on page 13 illustrates the schematic diagram of the test.




Figure 9 : Schematic diagram of the hot-generation filtration test

Methylene chloride is generated by hot evaporation using a STUART (Figure 10) and a HEIDOLPH peristaltic pump (Figure 11). The pump flow rate is checked with an OHAUS scale (Figure 12) connected to a computer.

The pump is set to a flow rate of 0.6 g/min to generate a concentration of 30 ppm. The equipment used during these tests is illustrated below:




Figure 10 : STUART Griddle (Model SB460, Serial Number R000101212)



Figure 11 : HEIDOLPH Peristaltic Pump (Pumpdrive Model 5101, No. 523-21010-00-2, Serial Number 111101371)




Figure 12 : OHAUS "Adventurer" Scale (AX201 model, serial number B709751985)

III. Presentation of the results

III. 1 Neutrodine Unisorb®


The breakthrough curves are presented on the Graph 3 and Graph 4.

According to the mass evaporated:



Graph 3 : Neutrodine Unisorb® Class 1 breakthrough curve as a function of evaporated mass

- According to the duration :



Graph 4 : Neutrodine Unisorb® Class 1 breakthrough curve as a function of time

According to these data, the 1% concentration of the occupational exposure limit value for methylene chloride according to the EPA (i.e. 2 ppm) reached downstream of the filtration system after the adsorption of 201.5 g of pollutant in 300 min.

These results are determined from the average of the two tests.

The Table 1 summarizes the results of the two tests.


	At 0.02 ppm	At 0.5 ppm
Evaporated mass (g)	170	222
	233	291
Average	201,5	256,5
Retention capacity (g/g)	0,0134	0,0175
	0,0183	0,0229
Average	0,0159	0,0202
Time (min)	270	350
	330	400
Average	300	375

Table 1 : Summary of the tests on Neutrodine Unisorb® in class 1

III.2 AS


The breakthrough curves are presented on the Graph 5 and Graph 6.

- According to the mass evaporated:



Graph 5 : AS class 1 breakthrough curve as a function of evaporated mass

- According to the duration :



Graph 6 : Class 1 AS breakthrough curve as a function of time

Based on these data, the 1% concentration of the occupational exposure limit value for methylene chloride according to the EPA (i.e. 2 ppm) is reached downstream of the filtration system after the adsorption of 123.3 g of pollutant in 146.5 min.

These results are determined from the average of the two tests.

The Table 2 summarizes the results of the two tests.


	At 0.02 ppm	At 0.5 ppm
Evaporated mass (g)	113,6	197
	133	301
Average	123,3	249
Retention capacity (g/g)	0,0118	0,0205
	0,0128	0,0291
Average	0,0123	0,0248
Time (min)	147	276
	146	328
Average	146,5	302

Table 2 : Summary of tests on AS in class 1

III.3 Neutrodine®


The breakthrough curves for both trials are presented on the Graph 7 and Graph 8.

According to the mass evaporated:



Graph 7 : Neutrodine® class 1 breakthrough curve as a function of evaporated mass

- According to the duration :



Graph 8 : Neutrodine® Class 1 breakthrough curve as a function of time

Based on these data, the 1% concentration of the occupational exposure limit value for methylene chloride according to the EPA (i.e. 2 ppm) is reached downstream of the filtration system after the adsorption of 56.7 g of pollutant in 71.5 min.

These results are determined from the average of the two tests.

The Table 3 summarizes the results of the two tests.

	At 0.02 ppm	At 0.5 ppm
Evaporated mass (g)	54	127
	59,4	145
Average	56,7	136
Retention capacity (g/g)	0,0036	0,0086
	0,0043	0,0105
Average	0,0040	0,0095
Time (min)	70	150
	73	180
Average	71,5	165

Table 3 : Summary of the tests on Neutrodine® in class 1

III.4 Summary

The following tables summarize the results of these filtration tests presented in this report.

At 0,02 ppm			
	2C Neutrodine Unisorb	2C AS	2C Neutrodine
Evaporated mass (g)	201,5	123,3	56,7
Time (min)	300	146,5	71,5

Table 4 : Summary of values at 0.02 ppm

At 0,5 ppm			
	2C Neutrodine Unisorb	2C AS	2C Neutrodine
Evaporated mass (g)	256,5	249	136
Time (min)	375	302	165

Table 5 : Summary of values at 0.5 ppm

Time at 0 ppm			
	2C Neutrodine Unisorb	2C AS	2C Neutrodine
Time (min)	274	128,5	54,5

Table 6 : Summary of times at 0 ppm

IV. Conclusion

Following new EPA regulations that restricted the TLV of methylene chloride to 2 ppm, new Class 1 filtration tests on Neutrodine Unisorb®, AS, and Neutrodine® were performed. TD/GC/MS achieved a detection limit < 1% US T: 0.02 ppm.

Thus, the retention rates have been updated on these different filtration columns:

- for Neutrodine Unisorb®, it is necessary to evaporate 201.5 g of pollutant in 300 min to reject 0.02 ppm.
- for AS, it is necessary to evaporate 123.3 g of pollutant in 146.5 min to reject 0.02 ppm.
- for Neutrodine®, it is necessary to evaporate 56.7 g of pollutant in 71.5 min to reject 0.02 ppm.

Neutrodine Unisorb® technology is the best suitable filtration technology to achieve the best lifetime before a release of 0.02 ppm methylene chloride.





About Erlab

The Erlab Research and Development laboratory

Since 1968, **Erlab** has been a specialist, inventor and world leader in **ductless, zero-emission filtering fume hoods for laboratories** to provide total safety in chemical handling.

1 Erlab filtration

We provide technologies to protect laboratory staff from inhaling chemicals. This is made possible thanks to our Research and Development (R&D) department, which has continuously improved our filtration technology for more than 50 years. That's why, in 2009, we invented the ERLAB ABOVE label for tried and tested filtration technology.

2 The AFNOR NF X 15-211: 2009 standard

Erlab's filtration technology conforms to the NF X 15-211: 2009 standard, the industry's most demanding standard for molecular filtration, developed by a committee of independent scientists and specialized manufacturers.

This text imposes performance criteria linked to:

- Filtration efficiency
- Containment efficiency
- Air face velocity
- Documentation: chemical listing

3 The ESP program

A set of three services included with the purchase of each device designed to ensure your safety.



eValQuest Risk analysis – Determination of protection needs – Determination of ergonomic needs.



ValiPass Certified installation – Total safety for handling.



ValiGuard Ongoing monitoring – Preventative and maintenance inspections – Device reconfiguration based on protection needs – Development of handling.

4 Flex technology

The combination of molecular and particulate filtration technologies allows a single device to meet laboratories' protection needs. This innovation from Erlab's R&D department offers unprecedented flexibility, versatility and value. A single device can be reconfigured over time and easily reassigned to other applications.

5 Smart technology

Smart technology is a simple and innovative means of communication that improves safety. This technology uses a light and sound signal to indicate the user's level of protection. The advantages of the technology are:

- 1/ Light pulsation: Real-time communication via LED light pulses intuitively alerts the user to the device's operating status.
- 2/ Simplicity: One-touch activation.
- 3/ Detection system: The exclusive detection system continuously monitors filtration performance.
- 4/ Built-in monitoring: This service provides direct access to the status, settings and history of your device.

France
+33 (0) 2 32 09 55 80 | ventes@erlab.net

Germany
0800 330 47 31 | export.north@erlab.net

United States
+1 800-964-4434 | info@erlab.com

United Kingdom
+44 (0) 1722 341 940 | export.north@erlab.net

China
+86 (0) 512 5781 4085 | sales.china@erlab.com.cn

Italy
+39 (0) 2 89 00 771 | export.south@erlab.net

Spain
+34 936 732 474 | export.south@erlab.net



erlab usa.erlab.com
You can breathe. iaq.erlab.com